



# E-xtreme

Make Essence to Extreme  
Smaller | Lighter | Quieter  
than ever

# Smaller

- 360° contra-angle rotation



- Miniature contra-angle for better visibility & easier posterior access.



Miniature contra-angle  
Super Mini head & thin neck




Normal head & neck

- Slim and compact design for maximum comfort and mobility





## Lighter



Less than 100 grams,  
reducing hand fatigue during treatment,  
providing the ultimate E-xtreme experience.



# Quieter

Eighteeth innovative motor,  
cutting down the noise from inside,  
allowing you to experience the tranquility.



Reciprocation function:  
five optional angles



Built-in file system



Model	: E-xtreme
Weight(grams)	: 99 g
Torque range	: 0.5 N.cm-4.0 N.cm
Power supply	: 100 V-240 V
Frequency	: 50 Hz / 60 Hz
Output	: DC 5 V/ 1 A